

## Multi-Stage Industrial Bioreactor Configuration



### Vessel

Capacity	1,000L-200,000L, customizable
Working Volume	30%-75%
Vessel Material	SUS316L/SUS304; smooth surface, no dead zones, easy to clean
Surface Finish	Inner surface mirror polished, $Ra \leq 0.4\mu\text{m}$ ; Outer surface mirror polished, $Ra \leq 0.6\mu\text{m}$
Bioreactor inlets	Vessel Lid: Standard port, charging port, feeding port, antifoam port, acid/alkali port, exhaust port, and sight glass port with LED light source Vessel Body: Sampling port; pH/DO and temperature ports; air inlet; discharge outlet; jacket water inlet and drain Liquid Level View: Long-strip liquid level sight glass inside the vessel
Air inlet	air inlet Mode: Deep aeration; air inlet pipe enters from the upper side wall of the vessel; internal air pipe is detachable Gas Distribution: Equipped with gas sparger for uniform gas

	<p>distribution without dead zones; easy to clean; meets high aerobic demand</p> <p>Backflow Protection: Equipped with non-return device to prevent material back-suction and contamination of inlet filter</p>
Discharge	The bioreactor is equipped with a bottom drain outlet featuring a sanitary, dead-zone-free design. It is easy to clean, supports repeated sterilization, and allows complete discharge of the culture broth from the vessel.
Sampling	A flush-mounted, dead-zone-free sterile sampling valve assembly is adopted, enabling multiple in-situ sterile sampling without compromising process sterility.

### Agitation System

Agitation Type	<p>The bioreactor is equipped with 3-stage, 6-blade impellers, with adjustable installation height. A mechanical defoaming impeller is also provided to effectively control foam formation during operation.</p> <p>The impeller type and configuration can be customized.</p>
----------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### Control System

Hardware Configuration	Siemens PLC + 15-inch LCD touch screen
Control Software	<p>Dedicated fermentation control software (independent intellectual property)</p> <p>Real-time and historical curves; supports 24-hour continuous monitoring and control</p> <p>Integrated fermentation process analysis software</p>
Automatic Sterilization	<p>Sterilization Mode: In-vessel steam sterilization; sterilization timing; over-temperature and over-pressure alarms</p> <p>Post-Sterilization Control: Automatic pressure holding and temperature control after sterilization to ensure sterility</p> <p>Filter Sterilization: Independent filter sterilization; set sterilization temperature and time; automatic drying after sterilization</p> <p>Sterilization Validation: Equipped with cold spot detection</p>

Temperature Control	<p>Display Accuracy: 0.3% FS±0.1°C  Control Accuracy: ±0.2°C  Sensor: WIKA Germany  Control Function : Dual-direction heating and cooling control; equipped with constant-temperature water tank  Data Management: Data logging, production curves, report analysis, remote data transmission</p>
pH	<p>Display Accuracy: 0.01 pH; measurement range: 0-12 pH  Sensor: METTLER (Germany) / HAMILTON  Control Function: Dual-direction acid and alkali addition; digital calibration  Data Management: Data logging, production curves, report analysis, remote data transmission</p>
DO	<p>Display Accuracy: 0.1%; range: 0-200%  Sensor: METTLER (Germany) / HAMILTON  Control Function: Cascade control; linked control with agitation speed and aeration  Data Management: Data logging, production curves, report analysis, remote data transmission</p>
Agitation Speed	<p>Speed Range: 0-300 rpm, customizable  Data Management: Data logging, production curves, report analysis, remote data transmission</p>
Gas Flow	<p>Control Function: Mass flow controller for automatic inlet gas control (maximum aeration up to 2.0 vvm; continuously adjustable according to process requirements)  Flow Meter: Vögtlin (Switzerland)/YOKOGAWA (Japan)</p>
Pressure	<p>Display Accuracy: 0.01 bar; range: 0-3 bar  Sensor: BD (Germany)  Control Function: Digital setpoint control for vessel pressure regulation  Data Management: Data logging, production curves, report analysis, remote data transmission</p>
Feeding	<p>Control Mode: Automatic control  Programmed Feeding: Time-based programmed feeding with adjustable variable speed  Batch Feeding: Set feed rate, single batch volume, and feeding cycle  Programmed Flow Feeding: Pre-set batch volume; multi-batch variable flow feeding based on time profile</p>